

Amendments to the Specification:

Please replace paragraph 3 on page 6 with the following amended paragraphs:

PCR Amplification and DNA sequencing. PCR amplification of the *DRD4* promoter polymorphism was conducted as described (18,19). The program OLIGO 6.0 was used to select primer pairs for the exon 1 polymorphism (20) (5'-TGGGCCCGCCGCATTCGT-3' (SEQ.ID.NO. 71) and 5'-GTGGGTGTATCGCCGAGGGA-3' (SEQ.ID.NO. 72); 661-nucleotide product) and the exon 3 VNTR (2) (5'-CGTACTGTGCGGCCTCAACGA-3' (SEQ.ID.NO. 73) and 5'-GACACAGCGCCTGCGTGATGT-3' (SEQ.ID.NO. 74); 705 nucleotide product for the 4R-allele). For some amplifications of the VNTR, primers described previously were used (2). The alternative primers were chosen farther from the VNTR, to minimize out-of-register hybridization during amplification. PCR reactions were conducted in 25 microliter volumes, containing 100ng genomic DNA, 200 micromolar dXTPs, 0.5 micromole of each primer, 1X PCR buffer (Qiagen), 1X Q-solution (Qiagen) and 0.625 units *Taq* DNA polymerase (Qiagen). Amplification was performed using Perkin-Elmer 9700 thermal cyclers. A 20 second, 96-degrees C hot start was used, followed by 40 cycles of 95 degrees C for 20 seconds and 68 degrees C for 1 minute. Following a 4-minute chase at 72-degrees C, excess primers were eliminated with 0.5 units of Shrimp Alkaline Phosphatase (SAP, Amersham Life Science), 0.1 unit of Exonuclease I (Exo I, Amersham Life Science) and 1X SAP buffer (Amersham Life Science). The SAP/Exo I reaction was carried out at 37 degrees C for 1 hour, followed by a 15-minute heat inactivation at 72-degrees C. The DNA from the SAP/Exo I reaction was used directly for DNA sequencing. For most individuals, the two allelic PCR products were first separated on 1.2-% agarose gels. DNA cycle sequencing was conducted by standard techniques, using ABI 377 and 3700 automated sequencers (21). DNA sequences of the *DRD4* haplotypes reported in this paper have been submitted to GenBank (Accession numbers AF395210 through AF395264).